Optional Autosampler for Convergys® X5 with integrated Barcode scanner and Mixer

Specifications

Sample capacity:	100 samples (10 individual
	racks for 10 x 10 sample
	tubes)
Autosampler	65 tests / hour
throughput:	
Tube Identification:	Via integrated barcode
	scanner, via on-screen or
	external keyboard (enter ID)
Sample preparation:	Via integrated sample mixer
Connection:	Plug and play connection via
	electronical connection port
Emergency sample:	Emergency break function
	allowing the manual
	measurement of emergency
	samples via Convergys X5
	sample rotor
Programming and	Easy-to-use, programming
Sampling:	and sampling menu via
	Convergys X5 interface
Operating conditions:	15°-30° C, 59-86° F
	(Optimal temperature is
	25 °C, 77 °F), 20%~ 80% RH
Dimensions:	$(W \times D \times H)$
	ca. 400x 300 x 210 mm
Net weight:	ca. 12 kg

Ordering Information

REF	Article Name	Description
1100-1700	Convergys® X5 Main Unit	Fully automatic 5-part WBC differential Hematology Analyzer
1100-1710	Convergys [®] X5 Autosampler	Autosampler for Convergys® X5 Hematology Analyzer
1100-1701	Convergys [®] Dil Diff D50 (20L)	Isotonic diluent
1100-1702	Convergys [®] Lyse-5P (5L)	Hemolysing Agent (WBC, HGB)
1100-1703	Convergys [®] Diff-5P (1L)	Hemolysing Agent (WBC, HGB, LYM, MON, NEU, EOS, BAS)
1100-1704	Convergys® Hypoclean CC (100ml)	Hypochloride cleaner
1100-1705	<i>Convergys®</i> check 5P Normal (3,0ml)	Blood control material, normal level
1100-1706	Convergys [®] check 5P High (3,0ml)	Blood control material, high level
1100-1707	Convergys [®] check 5P Low (3,0ml)	Blood control material, low level

* Convergent Technologies reserves the right to change any of the specifications without prior notice.

** Usage of original Convergent Technologies reagents is MANDATORY

*** Full specifications are available on request



Convergent Technologies GmbH & Co. KG Software Center 1, D-35037 Marburg, Germany Tel: +49 6421 8869948 | Fax: +49 6421 8869958 Homepage: www.convergent-technologies.de e-mail: info@convergent-technologies.de

CT_X5_VE

Convergys[®] X5 5-part WBC differential Hematology Analyzer



- 26 parameter 5-part WBC differential
- Laserlight scattering technology
- 100 µl sample volume of whole Blood
- 25 µl optional small sample module (SSM) for pediatric use
- 60 tests per hour throughput
- 65 tests per hour throughput with optional Autosampler
- Strict lysing process for best quality WBC differentiation
- Integrated clog prevention system
- 800x600 dots color graphic touch LCD
- Easy-to-use menu, based on Embedded Windows® XP® software
- Large data storage for 100.000 records





Convergys[®] X5 is one of the high-end models of Convergent-Technologies Hematology Analyzers range. The *Convergys*[®] X5 offers an optimal solution for hospitals, clinics and practices, which require a high throughput 26 parameter analyzer with *laser based optical measuring technology* for precise and accurate 5-part differential results for human samples. Our Hematology analyzers bring the convenience and accuracy of the reference laboratory right into your practice, offering maximum value at low running costs.


With the addition of the optional **small sample module (SSM*)** to the Convergys[®] X5, it is now possible to test samples having only a small blood volume in pediatric and neo-natal applications*.


***SSM** is validated for use with 3 tube types-BD Microtainer MAP, Sarstedt-Monovette and vacutainer to provide precise results on diluted samples. For use in open vial emergency mode only.

Convergys[®] X5 26 parameter 5-part WBC differential

High performance

- Complete 26-parameter CBC profile including the optical determination of the 5-part WBC differential count
- Autosampler with integrated bar code reader and mixer available
- Diagnostic flagging system
- Color touch-screen to monitor the results with large histograms and scattergrams
- Sophisticated operating menu: The instrument features an easy-to-use, logically constructed multilingual operating menu. It provides unlimited QC levels, L-J graphs and self-diagnostic functions to let you monitor reliability and accuracy
- The optical *laser diode* measurement technology ensures a longer lifetime and precise results with the special, high-tech flow cuvette
- *Convergys*[®]X5 ensures the data transfer with USB or RS232 port to ensure the link to your established LIMS (Laboratory Integrated Management System) for further processing of the measured data





Specifications of Convergys® X5

Measured Parameters:	CBC+5-part WBC Differential mode, 26 Para BAS%, EOS%, RBC, HGB, HCT, MCV, MCH,
Measuring Principles:	Volumetric Impedance Method for WBC ar Laser Light Scattering Technology for 5-p
Throughput:	60 tests / hour in the normal operations
Reagent System:	Isotonic Diluent, Hemolysing Agent (WBC BAS), Hypocleaner
Sampling Method:	Open and Closed tube System with Autom Primary Loops
Sample Volume:	100 μl of whole blood or 25 μl with small
Sample Types:	Human (general), Male, Female, Baby, Too
Tube Identification:	Over On-Screen or External Keyboard (Ent Over the Barcode Labels (Manual Barcode
Dilution Ratios:	WBC/BAS 1: 170, RBC/PLT 1: 21250, 4 DI
Chambers:	3 Chambers for Diluting Whole Blood and
Aperture Diameter:	70 μm (RBC/PLT), 80 μm (WBC/HGB)
HGB Measurement:	Integrated in WBC Chamber. Light source: Green LED with 540 nm Wav Detector: Light to Frequency Converter
TCU unit:	Transistors and Semiconductors are used 1 29° C for best quality WBC Differentiation
Optical Measurement:	Optical Measuring Head with closed Prote Light source: Semiconductor Laser Diode Sample Path: Quartz Flow Cell with Hydro Detector: Fibre-Optic Coupled PIN Si Phot
Auto-Alignment system:	Horizontal Calibration of Laser Beam Posi
Fine Calibration:	With Calibration Material (Polystyrene Mic
Clog Prevention:	High-Voltage Pulse on Aperture in Each A High-Pressure Back Flush of Aperture with
Cleaning Procedure:	High-Voltage Burst of the Aperture, High- Cleaner Reagent
Quality Control:	Unlimited QC Levels, QC Parameters inclue Parameters, 16- and 64-day Levey-Jennin
Calibration:	3, 5 or 7-Measurement SW Supported Aut
Flagging:	Warning Flags, Pathological (Diagnostic) Reagents Alert (Internal Buffers for Reage
User Interface:	800 x 600 Dots, Color Graphic Touch LCD, Embedded Windows® XP® Software
Multi-user Feature (advanced):	Multi-User Operation with Selective Privil
Languages Available:	English Menu with Support for German, S
Host Computer Interface:	USB Ports, Ethernet Port and RS-232 Seria
Data Back-up Method:	USB Mass Storage Device; External Lab Ma
Software Upgrade method:	USB Port, using USB Mass Storage Device
Data Storage capacity:	100.000 Records including Flagging, 4-Di
Data Processing:	Via C7 1.8 GHz Processor
Data Store:	Windows® XP® Embedded
Printer Interface:	Via USB Port to any Windows® XP® Compa
Display:	800 x 600 Dots, Color Graphic Touch LCD,
External Keyboard and Mouse:	Via PS/2 or USB Port
Barcode Reader:	Optional Manual Barcode Scanner via USB
Peripheral Ports:	USB 2.0 (4pc.), Ethernet, PS/2
Power Requirements:	100-120 or 200-240 VAC, 50-60Hz
Operating Conditions:	15° to 30° C, 59-86° F (Optimal tempera
Dimensions:	(W x D x H) ca. 515x 480 x 410 mm
Net weight:	ca. 35 kg

ameters: WBC, LYM, MON, NEU, BAS, EOS, LYM%, MON%, NEU% MCHC, RDWcv, RDWsd, PLT, P-LCC, P-LCR, PCT, MPV, PDWcv, PDWsd

nd RBC, Spectrophotometry for HGB, part WBC Differential

and 65 tests / hour with optional Autosamper

C, HGB), Hemolysing Agent (WBC, HGB, LYM, MON, NEU, EOS,

natic Sample Rotor, Ceramic Shear Valve with 3 Separated

Ill sample volume module (SSM) option

oddler and Child (Built-in Reference Ranges)

ter ID)

Scanner and/or Barcode Scanner in Autosampler)

IFF 1: 50

I Counting; 1:MIX, 1:RBC,1:WBC+HGB

velength

to maintain the Temperature Control Unit (TCU) Temperature at n and Maintain the Size of WBC Cells.

ective Housing

with 650 nm Wavelength and 7mW (Class IIIB)

o-Dynamic Focusing otodiodes with Internal Safety Interlock

sition.

icro-Particle or Polystyrene Microsphere, 5 μm)

Analysis Cycle; Chemical Cleaning;

th Cleaner Reagent

-Pressure Back Flush, Chemical Cleaning of the Aperture with

ude: Mean, \pm range, SD and CV for all Measured and Calculated ngs charts, Separate QC Database

comatic Mode and Manual (Factors) Calibration

Flags, Lab Limits (Normal Ranges) Jents), Instrument Alerts

Portrait Layout, Easy-to-Use, Menu Driven User Interface,

lege Levels, User Identification with ID and Password

Spanish, Croatian, Hungarian, Italian, Polish, Russian and Turkish ial Link.

anagement System, LIS via Ethernet or USB

)iff and BAS-Scattergrams, RBC and PLT Histograms

atible Printer

, Portrait Layout

3, Built-in Barcode Scanner in Optional Autosampler

ature is 25 °C, 77 °F), 20%~ 80% RH,









