Saturno 300 Plus

A true random access

CE Certified - ISO 9001 - IVD

CRONY®
INSTRUMENTS S.R.L.

Saturno 300 Plus

every patient is a stat

SATURNO 300 Plus is a fast random access analyzer for Clinical Chemistry and Immunoassay, where **EVERY PATIENT IS A STAT**

SATURNO 300 Plus is 24 hours on duty to run your daily routine, a single test, a profile, an electrolyte, an organ panel or a STAT with the same speed and efficiency.

State of the art technology with sophisticated user friendly software offers a new dimension in versatility, speed and reliability. It enables each user to personalize the analyzer to fit his particular needs and becomes an indispensable collaborator in any laboratory.

300 tests per hour plus 160 ISE electrolytes at the same time.

A total of 460 tests per hour

Operation Cycle

- Reagent is transferred into a reaction cuvette and brought to temperature prior to the addition of the sample.
- Cuvette with its reagent is immediately measured in the built-in photometer to obtain a Reagent/Cuvette Blank.
- The reaction cuvettes are immersed in a liquid bath at 37°C and are automatically washed for a non stop operation.
- Sample is transferred to ISE module, if required.
- Automatic print-out of patient report.

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
Il Programma

- F2** Will guide the operator step by step through all operations.
- F3** Can store thousands of patient data for future consultation on-line and to print patient reports.
- F4** To enter an unlimited number of user programmable methods.
- F5** To program an unlimited number of profiles or organ panels.
- F6** To program an unlimited number of standards, calibrators and multi-standard curves.
- F7** To enter an unlimited number of controls for the Quality Control program.
- F8** On-line control of all operating steps and procedures of the system. On-line trouble shooting.
- F9** To personalize the analyzer.

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Reagent Control

- 60 On-line reagents cooled or left at room temperature.
- On-line reagent inventory and display of the number of available tests for each reagent container.
- On-line reagent stability check.
- Bar code reader for reagent identification (optional).

Patienty entry



- 40 On-line primary tubes or sample cups with continuous sample loading and programming while analyzer is in operation.
- STATs will be processed within 6 minutes, without disrupting routine testing.
- Bi-directional interface for receiving and transmitting data to an external computer.
- Bar code reader for sample identification (optional).

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

Methodologies

- Open system, unlimited number of methods all user programmable.
- Only 300 µl of reagent per test.
- Automatic sample pre-dilution, if necessary.
- Two reagents chemistries and immunoassays.

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
60 days of Q.C. data displayed with actual values.

Typical Kinetic Measurement


Graph can be displayed for every sample.

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Saturno 300 Plus

Highlights

- Continuous patient programming while analyzer is in operation.
- The extractable reagent tray contains 60 on-line cooled reagents.
- Cost effective operation - no disposables.
- Walk away system - the reaction cuvettes are automatically washed and dried to assure a non stop, continuous operation.
- Sample and reagent dispensing probes with built-in liquid sensor and mixer.
- Check and control of reagent inventory and stability, flagged when reagent is missing or expired.
- Automatic repetition of critical or out of range results.
- I.S.E. module for simultaneous determination of Na, K and Cl (optional).
- Identification of sample and reagent with a built-in bar code reader (optional).
- Check of the quality of the washing and incubation water.
- Help on-line.

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Saturno 300Plus

A TRUE RANDOM ACCESS ANALYZER

OPERATING	True random access
ASSAY METHODS	Kinetic, End point, Fixed Time, IR-Kinetic, IR-End Point, IR Fixed Time, Differential
CALCULATION	Factor, Std, Point to Point, Linear, Quadric, Log-Logit 4, Log-Logit 5, Tetric, Lagrange, Sigmoid, Cubic, Cubic Spline, Exponential
REAGENTS ON LINE	60 Reagent on line
OPEN SYSTEM	Unlimited number of methods, standards, profiles, calibrators and controls user programmable
THROUGHPUT	More than 330 test/h + 160 electrolytes
SAMPLES TRAY	Extractable with 60 places, 40 are set aside for samples and 20 for Std, Ctri, Calibrators and STATs. Each position has a corresponding cup for automatic pre-dilution, if required.
SAMPLE VOLUME	From 3 to 70 μ l (increments of 1 μ l)
REAGENT VOLUME	R1 from 3 to 600 μ l (increments of 1 μ l) R2 from 3 to 597 μ l (increments of 1 μ l)
MINIMUM VOLUME	Only 300 μ l per test
CUVETTE VOLUME	Max 600 μ l
OPTICAL PATH	7 mm
PHOTOMETER	Multi-filter photometer with 8 interferential filters from 340 to 700 nm
ANALYTICAL PLATE	39 cylindrical cuvettes in quartz, immersed in a 37°C liquid bath. The built-in automatic washingsystem assures non-stop operation.
LAMP	Halogen 12V/20W - Pre-calibrated - Long life
DETECTOR	Silicon Photodiode
PHOT. LINEARITY	From -500 to 2,800 Abs \pm 1%
PHOT. SENSITIVITY	0,0005 Abs
ISE MODULE	For Na, K and Cl (optional)
BAR CODE READER	Identification of samples and reagents (optional)

COMPUTER

- Computer built-in - Pentium 1 GHz, 1 Gb RAM, HD >250Gb
- Windows XP Embedded
- Serial port
- USB port
- Bidirectional interface LIS
- Remote access for service (optional)

GENERAL

Power supply	230/115 Vac 50/60 Hz - 400W
Dimensions	103x67x45 cm
Weight	104 kg

Technical specifications are subject to change without prior notice

Via Tiburtina 661 - 00159 Rome (Italy)
Tel. +39 06 4396432 - +39 06 4381163 Fax +39 06 4390487
Web site: www.crony.com
E-mails: mail@crony.it - mail@crony.com - sales@crony.it

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