


# Gymna ShockMaster

## The radial shockwave therapy system

### Operational principle

A rapidly-moving projectile transfers its energy to the treatment head. The generated radial shockwaves are spread out and absorbed by the area of the body that is being treated.



### Effects

- Clear and immediate pain reduction
- Improvement of mobility (ROM)
- Removal of fibrosis and calcifications
- Stimulation of the metabolism

### Benefits

- The ideal therapy for chronic and complex musculoskeletal complaints
- Success rate of 80%
- Short treatment times: approx. 10 minutes
- 6 to 8 treatment sessions needed on average
- Realistic alternative to surgery
- Completes the range your practice can offer



# My best investments ever!

Gymna ShockMaster



Cryotherapy



Diathermy



Electrotherapy



Whole Body Vibration



Treatment Couches



Your dealer:

**gymna**<sup>®</sup>

[www.gymna.com](http://www.gymna.com)

GymnaUniphy NV Pasweg 6A | B-3740 Bilzen, België

Tel: +32 (0)89 510510 | Fax: +32 (0)89 210513

# Gymna ShockMaster 300

## Radial shockwave therapy anytime, anywhere.

# NEW!



**gymna**<sup>®</sup>

[www.shockmaster.eu](http://www.shockmaster.eu)

art nr. 38339 - ENG © Gymna Uniphy 2011



# Gymna ShockMaster 300

Radial shockwave therapy is an established and recognized form of treatment. The portable ShockMaster 300 means you are now able to carry out treatments anytime and anywhere. No restriction – just ‘plug and play’.



## Mobile

*Go on the road with this portable unit*

- Small and lightweight device: easy to install
- Built-in air compressor: no maintenance

## Comfort

*Ideal working conditions – the patient gets your undivided attention.*

- Perfect viewing angle for you and your patient
- Protective mat: safe and simple way of putting the hand piece down
- Silent air compressor: optimum communication between therapist and patient



D-ACTOR applicator



## ShockMaster 300

- Silent built-in compressor
- Max. frequency: 15 Hz
- Max. pressure: 4 bar
- Weight: 11 kg
- Dimensions (W x D x H): 42 x 38 x 18 cm
- Carrying case (optional)
- 15 mm treatment head
- D-Actor 20 mm treatment head
- 4 revisions
- 7" LCD full-colour touchscreen

## Time savings

*Save time by using the ShockMaster*

- DIY revision: saves time and money
- Colourful illuminated bar: clear representation of the treatment time
- Large operating buttons: changing parameters is quick and easy



Standard 15 mm applicator

Radial shockwaves

Tendinopathies



D-ACTOR applicator 20 mm

Radial shockwaves and vibrating pulses

Myofascial trigger points, muscle and connective tissue, fibrosis, deeper lying regions of muscle

see [www.shockmaster.eu](http://www.shockmaster.eu) for the other applicators

# Indications



**Shoulder pain, with or without calcification**  
Improvement of mobility.



**Plantar fasciopathy**  
Improved ability to take weight and reduction of pain.



**Myofascial trigger points**  
Disruption of permanently contracted sarcomeres in the muscle fibres.  
Reduction of radiating and local pain.

## Indications

- Tennis elbow/Lateral epicondylitis
- Golfer's elbow/Medial epicondylitis
- Dupuytren's contracture
- Ledderhose disease
- Heel spur/Plantar fasciopathy
- Shin splints/Anterior Tibial Syndrome
- Patellar tendinopathy
- Achilles tendon pain
- Myofascial pain syndrome
- and others...